

Aspirators

Aa



A763-010

Aspirator Bottles

- Tubular side arm for attachment of 12mm bore flexible tubing
- Unground neck to accept rubber stoppers

Catalogue No	Capacity, ml	Body diameter, mm	Maximum internal neck diameter, mm	Bottle height, mm
A763-006	250	73	18	131
A763-008	500	89	23	163
A763-010	1000	111	28	200
A763-012	2000	138	28	248

Bb


Cc

Dd

Ee

Ff

Gg



A763-042

Aspirator Bottles

- Complies with ISO 4796
- Ground side socket (size 29/32) to accept stopcock assemblies
- Unground top neck to accept rubber stoppers

Catalogue No	Capacity, ml	Body diameter, mm	Maximum internal neck diameter, mm	Bottle height, mm
A763-032	2000	138	28	248
A763-036	5000	181	44	318
A763-038	10000	227	59	398
A763-042	20000	288	59	492

Hh

Ii

Jj

Kk

Ll



A763-065

A763-063

Aspirator Stopcocks

For use with Aspirator bottles with a 29/32 side socket.

* Taken from cone shoulder, i.e. protrusion length.

Catalogue No	Description	Stopcock bore, mm	Overall length, mm
A763-063	Rotaflo® PTFE Key, 29/32	3	53
A763-065	Pyrex Glass Key, 29/32	4	95*

Mm

Nn

Oo

Pp

Qq

Rr

Ss

Tt

Uu

Vv

Ww

Xx

Yy

Zz



A763-086

A763-091

Aspirator Bottles, Large Capacity

- Ideal for storage of distilled water etc.
- Wide neck with simple glass cover
- Supplied with plastic connector and rubber/PTFE sealing ring
- Accepts stopcock A763-091 (listed below or glass tubing of 14mm outside diameter)

Catalogue No	Capacity, ml	Ext. neck diameter, mm	Total height, mm	Max diameter, mm
A763-080	5000	93	300	189
A763-082	10000	113	368	236

Aspirator Stopcock

Rotaflo tap with inert PTFE key- requires no lubrication.

Catalogue No	Stopcock bore, mm
A763-091	10

Aspirators, Narrow Neck, PP

- Octagonal shape with integral handle
- Polypropylene construction for high chemical resistance
- Autoclavable at 121°C

Catalogue No	Capacity, litres	Cap size, mm	Max. (h x w), mm
A763-512	5	45	270 x 220
A763-514	10	90	360 x 280
A763-516	20	90	420 x 350



Aspirators, Narrow Neck, HDPE

- Ideal for dispensing of common solvents, distilled water etc
- Fitted with carrying handles for ease of transport

Catalogue No	Capacity, litres	Neck i.d., mm
A763-531	5	49
A763-533	10	49
A763-535	25	64
A763-537	50	64



Aspirators, Wide Neck, HDPE

- Ideal for dispensing of common solvents, distilled water etc
- Easy fill wide aperture necks with water tight screw caps and carrying handles for ease of transport

Catalogue No	Capacity, litres	Neck i.d., mm	Height, mm
A763-568	30	240	415
A763-570	50	240	640
A763-572	100	240	770



Aspirators, LDPE

- Graduated to deliver in 5 litre increments
- Suitable for collecting distilled water, reagents and acids
- Wide shoulder handles

Catalogue No	Capacity, litres	Height, mm	Outer dia. mm	Closure size
A763-854	10	389	250	83B
A763-856	20	528	286	83B
A763-858	25	594	287	83B
A763-860	50	678	379	83B



Aa

Bb

Cc

Dd

Ee

Ff

Gg

Hh

Ii

Jj

Kk

Ll

Mm

Nn

Oo

Pp

Qq

Rr

Ss

Tt

Uu

Vv

Ww

Xx

Yy

Zz