

Basins

Aa

Bb

Cc

Dd

Ee

Ff

Gg

Hh

Ii

Jj

Kk

Ll

Mm

Nn

Oo

Pp

Qq

Rr

Ss

Tt

Uu

Vv


Ww

Xx

Yy

Zz

B073-124



Basins, Evaporating, Round Bottom, with Spout

Borosilicate glass, Simax.

Catalogue No	Capacity, ml	Diameter, mm	Height, mm	Pack qty.
B073-120	50	64	29	10
B073-124	100	85	39	10
B073-128	250	105	49	10

B073-135



Basins, Evaporating, Flat Bottom, with Spout

Borosilicate glass, Simax.

Catalogue No	Capacity, ml	Diameter, mm	Height, mm	Pack qty.
B073-131	45	60	30	10
B073-135	90	80	35	10
B073-139	170	95	40	10

Basins, Evaporating

Pyrex®

- Manufactured from Pyrex® borosilicate glass
- Suitable for use up to 400°C
- Chemically and thermally resistant

Catalogue No	Capacity, ml	Diameter, mm	Height, mm
B073-213	100	70	40
B073-215	150	80	45
B073-217	300	95	55
B073-219	500	115	65
B073-221	900	140	75
B073-223	2000	190	90

B073-113



Basins, Evaporating, with Spout

Pyrex®

- Manufactured from Pyrex® borosilicate glass
- High resistance to chemical attack and thermal shock
- Flat base curved sides with spout

Catalogue No	Capacity, ml	Diameter, mm	Height, mm
B073-242	15	50	25
B073-244	45	60	30
B073-246	90	80	45
B073-248	170	95	55
B073-250	320	115	65
B073-252	600	140	80

B073-242



Basins, Evaporating, Porcelain, Round Bottom, with Spout

Catalogue No	Diameter, mm	Capacity, ml	Height, mm	Pack qty.
B073-660	65	36	28	10
B073-662	80	72	33	10
B073-664	85	107	38	10
B073-666	95	140	42	10
B073-668	105	190	45	10
B073-670	115	232	47	10
B073-672	130	390	55	10
B073-674	150	590	66	10

B073-662

