

Brushes for Test Tubes

With fan-shaped end, on galvanised steel wire handle.

Catalogue No	For tube diameter, mm	Brush head dia, mm	Pack qty.
Bristle			
B677-410	10-12	11	10
B677-416	16-25	29	10
Nylon			
B677-422	10-12	11	10
B677-425	16-25	29	10

Brushes for Test Tubes

With cotton yarn tip on galvanised steel wire handle.

Catalogue No	For tube diameter, mm	Brush head dia, mm	Pack qty.
Bristle			
B677-442	12-16	24-13 (tapered)	10
B677-446	16-25	35-25 (tapered)	10
Nylon			
B677-454	12-16	24-13 (tapered)	10
B677-458	16-25	35-25 (tapered)	10

Brushes for Beakers

Bristle with 40cm wooden handle.

Catalogue No	Bristle dimensions (wxdx), mm	Pack qty.
B677-462	65 x 45 x 150	10

Brushes for Burettes

White bristle on galvanised steel wire handle.

Catalogue No	For burette capacity, ml	Brush head dia, mm	Length, mm	Pack qty.
Bristle				
B677-512	25	16	610	10
B677-514	50	19	760	10
B677-518	100	25	915	10
Nylon				
B677-524	25	16	610	10
B677-526	50	19	760	10
B677-528	100	25	915	10

Brushes for Pipettes and Tubes


Head size 125 x 10 mm. Overall length 700 mm. For pipettes and tubes from 8-12 mm diameter.

Catalogue No	Brush type	Pack qty.
B677-551	Bristle	10
B677-553	Nylon	10

Brushes General Purpose

Suitable for use with measuring cylinders, flasks, bottles and beakers diameter from 50-100mm. With fan-shaped end on galvanised wire handle.

Catalogue No	Brush type
B677-572	Bristle
B677-574	Nylon



Aa

Bb

Cc

Dd

Ee

Ff

Gg

Hh

Ii

Jj

Kk

Ll

Mm

Nn

Oo

Pp

Qq

Rr

Ss

Tt

Uu

Vv

Ww

Xx

Yy

Zz