

# Burners

Aa

Bb

Cc

Dd

Ee

Ff

Gg

Hh

Ii

Jj

Kk

Ll

Mm

Nn

Oo

Pp

Qq

Rr

Ss

Tt

Uu

Vv

Ww

Xx

Yy

Zz

B757-100



## Bunsen Burner, 13mm Outer Diameter

All burners have a heavy pressure diecast base, 80mm diameter, painted in a blue chemically resistant acrylic finish.

Catalogue No	Overall height, mm	For gas type
<b>B757-100</b>	145	LPG
<b>B757-104</b>	145	Natural

B757-122



## Bunsen Burner, Control Stopcock, Pilot Ignition Flame

Catalogue No	Overall height, mm	For gas type
<b>B757-122</b>	171	LPG
<b>B757-126</b>	171	Natural

B757-141



## Bunsen Burner, Meker Pattern

Manufactured from Nickel plated Brass. The burner head is 24mm diameter, and fitted with a perforated steel head, which gives off multiple small but very hot flames.

Catalogue No	Overall height, mm	For gas type
<b>B757-141</b>	162	LPG
<b>B757-143</b>	162	Natural

B757-160



## Burner, Universal

Useful for sterilising laboratory instruments, the burner emits an instant safety flame, with a push button ignition which incorporates a LOCK/ON mechanism. The unit switches off when the users hand is removed. Simple gas fill using butane. Long flint life. Powered by a 1.5 V battery.

Catalogue No	<b>B757-160</b>
--------------	-----------------

B757-176



## Cartridge Burner

Portable burner for use with a 220gm Butane/Propane mix gas cartridge.

- Typical burning time of the 220gm cartridge is 2.5 hours
- Overall height when mounted on cartridge is 190mm

Catalogue No	<b>B757-176</b>
--------------	-----------------

B757-184



## Spirit/Alcohol Lamps

- For use with methylated spirit
- Glass reservoir, 100ml, complete with wick and lid

Catalogue No	<b>B757-184</b>
--------------	-----------------

## Accessories: Tubing for Bunsen Burners

Flame retardant neoprene tubing, with reinforced ends.

Catalogue No	Length, mm
<b>B757-190</b>	450
<b>B757-192</b>	600
<b>B757-194</b>	900