

Aa
Bb
Cc
Dd
Ee
Ff
Gg
Hh
Ii
Jj
Kk
Ll
Mm
Nn
Oo
Pp
Qq
Rr
Ss
Tt
Uu
Vv
Ww
Xx
Yy
Zz

C257-100



C257-112

Centrifuge, Micro, Haematocrit, Economy Model



Brushless drive centrifuge.

- Variable speed
- Digital control of parameters
- 1-90 seconds pulse (quick spin)
- Speed, lid open & imbalance detectors
- CE Marked

Catalogue No	Max Speed, rpm	Max. rcf	Dims. (hwxwd), cm	Weight, Kg
C257-100	12000	13500g	23 x 26 x 39	14

Accessory Rotors

Catalogue No	Part No	Description	Max. speed, rpm
C257-112	5401	Rotor for 24 x Capillary Tubes, 75mm	12000
C257-114	5404	Rotor for 24 x 1.5-2ml Tubes	12000

C257-120



C257-127

Centrifuge, Micro, Variable Speed, Economy Model



Brushless drive centrifuge.

- Variable speed
- Variable time, pulse and hold
- Digital control of parameters
- Speed, lid open & imbalance detectors
- CE Marked

Catalogue No	Speed, rpm	Max rcf	Dims. (hwxwd), cm	Weight, Kg
C257-120	6000	3800g	23 x 26 x 39	14

Accessory Rotors

Catalogue No	Part No	Description	Max. speed, rpm
C257-125	5501	Swing out rotor for 8 x 4/5ml tubes	3750
C257-127	5420	Rotor for 12 x 15ml Tubes	6000
C257-129	5435	Rotor for 6 x 50ml Tubes	6000

C257-140



C257-154

Centrifuge, Variable Speed, Maximum Capacity 8 x 50ml, Economy Model



Brushless drive centrifuge.

- Variable speed
- Variable time, pulse and hold
- Digital control of parameters
- Speed, lid open & imbalance detectors
- CE Marked

Catalogue No	Speed, rpm	Max. rcf	Dims. (hwxwd), cm	Weight, Kg
C257-140	6000	4500g	26 x 33 x 41	17

Accessory Rotors

Catalogue No	Part No	Description	Max. speed, rpm
C257-151	5431	Rotor for 24 x 15ml Tubes	6000
C257-154	5434	Rotor for 8 x 50ml Tubes	6000

Requires adapters, please turn to page 64