

Filter Discs

Aa

Bb

Cc

Dd

Ee

Ff

Gg

Hh

Ii

Jj

Kk

Ll

Mm

Nn

Oo

Pp

Qq

Rr

Ss

Tt

Uu

Vv

Ww

Xx

Yy

Zz



Filter Discs

- Complies with ISO 4793
- For use by glassblowers constructing special apparatus
- Manufactured from Pyrex borosilicate glass for resistance to heat and chemical corrosion
- Close control of pore size
- Accurate and consistent disc dimensions
- High mechanical strength to withstand full vacuum
- Convenient packaging

Catalogue No	Porosity grade	Nominal diameter, mm	Pack qty.
F010-010	0	10	10
F010-020	0	20	10
F010-030	0	30	10
F010-040	0	40	10
F010-050	0	50	10
F010-065	0	65	1
F010-090	0	90	1
F011-010	1	10	10
F011-020	1	20	10
F011-030	1	30	10
F011-040	1	40	10
F011-050	1	50	10
F011-065	1	65	1
F011-090	1	90	1
F011-120	1	120	1
F012-010	2	10	10
F012-020	2	20	10
F012-030	2	30	10
F012-040	2	40	10
F012-050	2	50	10
F012-065	2	65	1
F012-090	2	90	1
F012-120	2	120	1
F013-010	3	10	10
F013-020	3	20	10
F013-030	3	30	10
F013-040	3	40	10
F013-050	3	50	10
F013-065	3	65	1
F013-090	3	90	1
F013-120	3	120	1
F014-010	4	10	10
F014-020	4	20	10
F014-030	4	30	10
F014-040	4	40	10
F014-050	4	50	10
F014-065	4	65	1
F014-090	4	90	1
F014-120	4	120	1