

Aa
Bb
Cc
Dd
Ee
Ff
Gg
Hh
Ii
Jj
Kk
Ll
Mm
Nn
Oo
Pp
Qq
Rr
Ss
Tt
Uu
Vv
Ww
Xx
Yy
Zz



Laboratory Refrigerators, Freezers and Fridge/Freezers



Economically priced range of refrigerators and freezers designed for the safe storage of samples and chemicals.

- Temperature display panel
- Spark free interior reduce the risk of internal explosion
- Flexible door storage with adjustable trays, lockable for product security
- CFC and HFC free refrigeration system and insulation
- Reversible door on all upright models
- Freezers are suitable for holding insulin, serum & plasma at correct temperatures

R226-023



Refrigerator Specifications

Catalogue No	R226-023	R226-026	R226-028
Model	LR207	LR507	LR907
Capacity, litres	67	145	285
Defrost		Automatic	
Temperature range, °C		+2 to +8	
Shelves/Drawers	2/1	3/1	5/1
External dimensions (wxdxh), mm	502 x 540 x 660	558 x 600 x 863	595 x 600 x 1550
Weight, Kg	27	34	81.5
Electrical supply		230V / 50Hz	

R226-045



Fridge/Freezer Specifications

Catalogue No	R226-042	R226-045	R226-049
Model	ISR27	ISR57	IST47
Capacity, litres (Fridge)	75	110	157
Capacity, litres (Freezer)	5.5	17	89
Fridge Defrost		Automatic	
Freezer Defrost		Manual	
Temperature range, °C (Fridge)		+2 to +8	
Temperature range, °C (Freezer)		To -18	
Shelves/Drawers (Fridge)	2/1	2/1	3/1
Shelves/Drawers (Freezer)	0	0	3
External dimensions (wxdxh), mm	500 x 515 x 636	558 x 600 x 863	560 x 600 x 1572
Weight, Kg	27	36	54
Electrical supply		230v / 50Hz	

R226-088



Freezer Specifications

Catalogue No	R226-082	R226-085	R226-088
Model	ISU27	ISU57	ISU97
Capacity, litres	38	94	202
Defrost		Manual	
Temperature range, °C		-18	
Shelves/Drawers	3	4	6
External dimensions (wxdxh), mm	500 x 515 x 636	558 x 600 x 863	595 x 600 x 1550
Weight, Kg	17	38	75
Electrical supply		230v 50Hz	