


## Drum, Bottle Holder, STR4/4

stuart

S308-538



Consists of a four segment Perspex® cradle. Each segment is fitted with a fully adjustable Velcro® strap to hold bottles and other containers up to 120mm diameter and up to 300mm in length.



Catalogue No	Description
S308-538	Drum with four segment cradles holding up to four bottles.

## Sieves

International Standard Test sieves, ISO 3310/1. All 200mm wire mesh sieves have a stainless steel frame and mesh. The 200mm perforated plate sieves also have a stainless steel frame. The sieves have a depth of 50mm.

Catalogue No	Aperture Size	Catalogue No	Aperture Size	Catalogue No	Aperture Size
S340-020	20 µm	S340-224	224 µm	S341-200	2.00 mm
S340-025	25 µm	S340-250	250 µm	S341-224	2.24 mm
S340-032	32 µm	S340-280	280 µm	S341-236	2.36 mm
S340-036	36 µm	S340-300	300 µm	S341-250	2.50 mm
S340-038	38 µm	S340-315	315 µm	S341-280	2.80 mm
S340-040	40 µm	S340-355	355 µm	S341-315	3.15 mm
S340-045	45 µm	S340-400	400 µm	S341-335	3.35 mm
S340-050	50 µm	S340-425	425 µm	S341-355	3.55 mm
S340-053	53 µm	S340-450	450 µm	S341-400	4.00 mm
S340-056	56 µm	S340-500	500 µm	S341-450	4.50 mm
S340-063	63 µm	S340-560	560 µm	S341-475	4.75 mm
S340-071	71 µm	S340-600	600 µm	S341-500	5.00 mm
S340-075	75 µm	S340-630	630 µm	S341-560	5.60 mm
S340-080	80 µm	S340-710	710 µm	S341-630	6.30 mm
S340-090	90 µm	S340-800	800 µm	S341-670	6.70 mm
S340-100	100 µm	S340-850	850 µm	S341-710	7.10 mm
S340-106	106 µm	S340-900	900 µm	S341-800	8.00 mm
S340-112	112 µm	S341-100	1.00 mm	S341-900	9.00 mm
S340-125	125 µm	S341-112	1.12 mm	S341-950	9.50 mm
S340-140	140 µm	S341-118	1.18 mm	S342-100	10.00 mm
S340-150	150 µm	S341-125	1.25 mm	S342-112	11.20 mm
S340-160	160 µm	S341-140	1.40 mm	S342-125	12.50 mm
S340-180	180 µm	S341-160	1.60 mm	S342-132	13.20 mm
S340-200	200 µm	S341-170	1.70 mm	S342-140	14.00 mm
S340-212	212 µm	S341-180	1.80 mm	S342-160	16.00 mm



## Accessories

Catalogue No	Description
S342-220	Lid
S342-224	Receiver
S342-228	Receiver/separator

## Sieve Shaker, A060

This model is activated by an electric motor.

- Holds up to 8 sieves dia. 200 mm or 8" plus pan and lid

Catalogue No	S342-330
Model	A060, Sieve Shaker
Power supply	220-240V, 50Hz, 1ph, 110W
Dimensions, mm	320 x 380 x 800
Weight, Kg	20

S342-330



Aa

Bb

Cc

Dd

Ee

Ff

Gg

Hh

Ii

Jj

Kk

Ll

Mm

Nn

Oo

Pp

Qq

Rr

**Ss**

Tt

Uu

Vv

Ww

Xx

Yy

Zz