

Aa
Bb
Cc
Dd
Ee
Ff
Gg
Hh
Ii
Jj
Kk
Ll
Mm
Nn
Oo
Pp
Qq
Rr
Ss
Tt
Uu
Vv
Ww
Xx
Yy
Zz



Bench Protection, Benchguard®

A highly absorbent white paper, which is plastic coated on the reverse side to protect work surfaces against accidental spillages. The top layer instantly absorbs the spillage, whilst the bottom layer prevents any seepage through to the bench. Benchguard® Extra has double the absorbency for use in more demanding applications.

Catalogue No	Description	Size	Absorption, ml/m ²	Pack qty.
S078-484	Benchguard®	50m x 49cm	400	Roll
S078-486	Benchguard®	60cm x 49cm	400	50 sheets
S078-488	Benchguard®	50m x 92cm	400	Roll
S078-492	Benchguard® Extra	50m x 49cm	800	Roll
S078-496	Benchguard® Extra	60cm x 49cm	800	50 Sheets



Safety Shields, Hinged

Transparent polycarbonate shield, supplied as three panels, joined by a polypropylene hinge. The shields are supported by a pair of heavy plastic coated metal feet, which can be clamped or bolted to the bench for greater stability.

Catalogue No	Height, mm	Outside panel width, mm	Centre panel width, mm	Total width, mm
S083-002	380	180	230	590
S083-004	610	180	230	590
S083-006	610	226	300	750



Safety Shield

Manufactured from a single piece of 4mm transparent polycarbonate. The shield is supported by a pair of heavy plastic coated metal feet, which can be clamped or bolted to the bench for greater stability.

Catalogue No	Height, mm	Outside panel width, mm	Centre panel width, mm
S083-102	608	178	230

Sound Level, Auto Ranging

Sound level meter with an electret condenser microphone giving virtually instant readings updated every 0.5 seconds on a large LCD display. Supplied with battery.

- 4 measuring ranges, 40-70dB, 60-90dB, 80-110dB, 100-130dB
- Auto ranging or manual range selection
- Quasi analogue bargraph, 1dB steps updated every 40ms
- Min/max memory since power on or last reset
- Auto power off

S085-510



Catalogue No	S085-510
Range	40 to 130dB ('A' weighting), 45 to 130dB ('C' weighting)
Resolution	0.1dB display, 1dB bargraph
Accuracy	±2dB
Power	9V PP3 battery
Dimensions (l x w x d), mm	180 x 72 x 30
Weight, g	184

Fosse Liquitrol: Spillage Binder, Uni-Safe

An odourless light green powder, which is capable of binding and absorbing all liquid spills, oil and chemical (except Hydrofluoric Acid). The powder changes colour to yellow when in contact with Acids, and green with Alkalis. It binds most liquids into a non-flowing gel and leaves the surface smear free, and

S092-335



Catalogue No	Pack Size, kg	Optimum Sorbent Capacity
S092-331	1	25 Litres Water, 12 Litres Oil
S092-335	5	125 Litres Water, 60 Litres Oil