

# Sieving, Spatulas

Aa  
Bb  
Cc  
Dd  
Ee  
Ff  
Gg  
Hh  
Ii  
Jj  
Kk  
Ll  
Mm  
Nn  
Oo  
Pp  
Qq  
Rr  
**Ss**  
Tt  
Uu  
Vv  
Ww  
Xx  
Yy  
Zz




## Sieve Shaker, A059-02

This model is activated by electromagnetic impulses and has a triple vibrating action (vertical, lateral and rotational). For use with sieves dia. 200, 250, 300 and 315mm, 8"-12". It is recommended for accurate sieving tests and also on fine materials.

- Holds up to 10 sieves
- Digital control panel
- Sieving time from 1-999 minutes
- Pause between one vibration and the following one

Catalogue No	<b>S342-335</b>
Model	A059-02 Electromagnetic sieve shaker
Power supply	220-240V, 50/60Hz, 1ph, 450/750W
Dimensions, mm	380 x 440 x 1080
Weight, Kg	65



## Spatulas, Polypropylene

- Good chemical and heat resistance
- Autoclavable

Catalogue No	Length, mm	Pack qty.
<b>S365-060</b>	100	10
<b>S365-062</b>	152	10
<b>S365-064</b>	200	10
<b>S365-066</b>	252	10



## Spatulas, Stainless Steel, Palette Knife

Parallel sided blade, with black polypropylene handle, autoclavable.

Catalogue No	Blade length, mm
<b>S365-211</b>	105
<b>S365-213</b>	150
<b>S365-215</b>	200
<b>S365-217</b>	250
<b>S365-219</b>	300



## Spatula, Chattaway, Nickel

Catalogue No	Length, mm	Pack qty.
<b>S365-232</b>	100	5
<b>S365-234</b>	150	5
<b>S365-236</b>	200	5

## Spatula, Micro Chattaway, Nickel

Catalogue No	Length, mm	Pack qty.
<b>S365-265</b>	100	5
<b>S365-267</b>	150	5
<b>S365-269</b>	200	5

## Spatula, Chattaway, Stainless Steel

Catalogue No	Length, mm	Pack qty.
<b>S365-310</b>	100	5
<b>S365-312</b>	150	5
<b>S365-314</b>	200	5

## Spatula, Micro Chattaway, Stainless Steel

Catalogue No	Length, mm	Pack qty.
<b>S365-344</b>	100	5
<b>S365-346</b>	150	5
<b>S365-348</b>	200	5

## Spatula, Spoon End, Stainless Steel

One end spoon, and the other end flat.

Catalogue No	Length, mm	Pack qty.
<b>S365-361</b>	120	10
<b>S365-366</b>	150	10
<b>S365-369</b>	210	10

## Spatula, Micro Spoon End, Stainless Steel

One end spoon, and the other end flat.

Catalogue No	Length, mm	Pack qty.
<b>S365-374</b>	100	10
<b>S365-378</b>	150	10

## Spatula, Trulla, Stainless Steel

Trulla Pattern, 19mm wide at one end tapering to a 10mm scoop at the other.

Catalogue No	Length, mm	Pack qty.
<b>S365-386</b>	175	5

## Spatula, Nuffield Pattern, Stainless Steel

Catalogue No	Length, mm	Pack qty.
<b>S365-392</b>	140	10

## Spatulas, Glass Fibre, Nylon Filled

A range of robust, general purpose spatulas with excellent chemical resistance.

Catalogue No	Type	Length, mm	Capacity of spoon, ml	Pack qty.
<b>S365-464</b>	Flat ends	150	-	10
<b>S365-466</b>	Flat ends	180	-	10
<b>S365-472</b>	Spoon flat	180	1.5	10
<b>S365-474</b>	Spoon flat	210	1.8	10



## Note:

We supply an extensive range of equipment that may not feature in this catalogue. If you can't find a particular item please call for specifications and pricing information.

Aa

Bb

Cc

Dd

Ee

Ff

Gg

Hh

Ii

Jj

Kk

Ll

Mm

Nn

Oo

Pp

Qq

Rr

**Ss**

Tt

Uu

Vv

Ww

Xx

Yy

Zz