


Burette Holder

The burette is held in a three point grip between two formed stainless steel stops and an eccentric cam. The holder is fitted with a universal bosshead that may be clamped to either vertical or horizontal rods.

Catalogue No	Type
S413-165	Single
S413-168	Double



Retort Rings, Plain Stem

Zinc plated mild steel. Plain 8mm diameter stem, 140mm long.

Catalogue No	Internal diameter, mm
S413-240	56
S413-244	100
S413-248	150



Retort Rings with Bosshead

Zinc plated mild steel. The stem length to the centre of the bosshead is approximately 60mm.

Catalogue No	Internal diameter, mm
S413-251	56
S413-255	100
S413-259	150

Stands, Burette

Polished hardwood stand with plastic clamping screws. Rods are 460 x 16mm approximately.

Catalogue No	Type	Base size, mm
S420-264	Single	200 x 100
S420-268	Double	290 x 110



Laboratory Jacks

Stainless steel operating screws and spindles, with plates manufactured from an aluminium extrusion. The three larger sizes, all have a 10mm tapped hole in the right rear corner to take a standard retort stand rod.

Catalogue No	Plate dimensions, mm	Max height, mm	Max load, kg
S420-532	100 x 100	140	30
S420-534	150 x 150	275	60
S420-536	200 x 200	275	80
S420-538	320 x 260	470	90



Tripod Stand, Triangular, Stainless Steel

Catalogue No	Length of side, mm	Height, mm
S420-743	120	150
S420-745	120	210
S420-747	150	210



Aa

Bb

Cc

Dd

Ee

Ff

Gg

Hh

Ii

Jj

Kk

Ll

Mm

Nn

Oo

Pp

Qq

Rr

Ss

Tt

Uu

Vv

Ww

Xx

Yy

Zz

Aa

Bb

Cc

Dd

Ee

Ff

Gg

Hh

Ii

Jj

Kk

Ll

Mm

Nn

Oo

Pp

Qq

Rr

Ss

Tt

Uu

Vv

Ww

Xx

Yy

Zz

S442-154




Rack, Test Tube, Aluminium

Stackable Z-pattern with eight 13mm, six 19mm and five 28mm holes.

Catalogue No	Description
S442-154	Z Rack

Test Tube Racks

- High-grade stainless steel construction
- Dimensions of rack: 270mm long x 70mm wide, height to top of lifting handle: 138mm



Catalogue No	Description
B139-313	Test tube rack, stainless steel, 36 holes x 13mm dia
B139-317	Test tube rack, stainless steel, 26 holes x 17mm dia
B139-319	Test tube rack, stainless steel, 18 holes x 19mm dia
B139-326	Test tube rack, stainless steel, 16 holes x 26mm dia
B139-332	Test tube rack, stainless steel, 12 holes x 32mm dia
<i>Max number of racks per bath</i>	
Bath Size	4 litre=1 rack 8 litre=2 racks 14 litre=4 racks 22/28 litre=6 racks

Rack, Microcentrifuge Tube, Polypropylene

- Supplied in self assembly flat packs
- Will sink when used in a water bath, and can be stacked when not in use
- Autoclavable to 121°C
- Supplied in packs of 5

Catalogue No	Colour	Size (l x w x h), mm	No. places	For tube capacity, ml
S442-212	White	265 x 126 x 38	128	1.5
S442-214	Blue	265 x 126 x 38	128	1.5



Rack, Test Tube, Kit Form, Polypropylene

- Test tube racks supplied in self assembly flat packs
- Autoclavable to 121°C and will sink in water baths
- Alpha numeric place coding
- Supplied in packs of 5

Catalogue No	Colour	Size (l x w x h), mm	No. places	Max. tube dia. mm
S442-241	White	265 x 126 x 75	84	13
S442-243	Blue	265 x 126 x 75	84	13
S442-252	White	265 x 126 x 75	55	18
S442-254	Blue	265 x 126 x 75	55	18
S442-263	White	265 x 126 x 75	40	21
S442-265	Blue	265 x 126 x 75	40	21
S442-270	White	265 x 126 x 88	32	26
S442-272	Blue	265 x 126 x 88	32	26
S442-281	White	265 x 126 x 88	21	31
S442-283	Blue	265 x 126 x 88	21	31

Test Tube Racks, Unwire™

- Made of tough chemical resistant ResMer™
- Range of colours available
- Autoclavable
- Stackable when empty

Catalogue No	Tube size, mm	Size (lwxh), mm	Colour	Array
S442-521	13	200 x 102 x 57	White	6 x 12
S442-523	16	248 x 127 x 70	White	6 x 12
S442-525	20	250 x 102 x 83	White	4 x 10
S442-527	25	300 x 121 x 92	White	4 x 10
S442-529	30	283 x 108 x 83	White	3 x 8
S442-530	13	200 x 102 x 57	Orange	6 x 12
S442-532	16	248 x 127 x 70	Orange	6 x 12
S442-534	20	250 x 102 x 83	Orange	4 x 10
S442-536	25	300 x 121 x 92	Orange	4 x 10
S442-538	30	283 x 108 x 83	Orange	3 x 8
S442-541	13	200 x 102 x 57	Yellow	6 x 12
S442-543	16	248 x 127 x 70	Yellow	6 x 12
S442-545	20	250 x 102 x 83	Yellow	4 x 10
S442-547	25	300 x 121 x 92	Yellow	4 x 10
S442-549	30	283 x 108 x 83	Yellow	3 x 8
S442-550	13	200 x 102 x 57	Blue	6 x 12
S442-552	16	248 x 127 x 70	Blue	6 x 12
S442-554	20	250 x 102 x 83	Blue	4 x 10
S442-556	25	300 x 121 x 92	Blue	4 x 10
S442-558	30	283 x 108 x 83	Blue	3 x 8
S442-561	13	200 x 102 x 57	Green	6 x 12
S442-563	16	248 x 127 x 70	Green	6 x 12
S442-565	20	250 x 102 x 83	Green	4 x 10
S442-567	25	300 x 121 x 92	Green	4 x 10
S442-569	30	283 x 108 x 83	Green	3 x 8
S442-570	13	200 x 102 x 57	Red	6 x 12
S442-572	16	248 x 127 x 70	Red	6 x 12
S442-574	20	250 x 102 x 83	Red	4 x 10
S442-576	25	300 x 121 x 92	Red	4 x 10
S442-578	30	283 x 108 x 83	Red	3 x 8



...the affordable solution

- Aa
- Bb
- Cc
- Dd
- Ee
- Ff
- Gg
- Hh
- Ii
- Jj
- Kk
- Ll
- Mm
- Nn
- Oo
- Pp
- Qq
- Rr
- Ss**
- Tt
- Uu
- Vv
- Ww
- Xx
- Yy
- Zz