

Aa  
Bb  
Cc  
Dd  
Ee  
Ff  
Gg  
Hh  
Ii  
Jj  
Kk  
Ll  
Mm  
Nn  
Oo  
Pp  
Qq  
Rr  
**Ss**  
Tt  
Uu  
Vv  
Ww  
Xx  
Yy  
Zz

S512-270



## Stirrers, Magnetic, Portable, SM27

stuart

Portable magnetic stirrer for use in glove boxes, incubators or in the field.

- Housed in chemically resistant ABS case
- Battery powered or mains operated
- Supplied complete with 4 x LR20 batteries and a 25mm PTFE stirrer bar
- With BioCote® antimicrobial protection

Catalogue No	<b>S512-270</b>
Capacity, litres	1.5
Maximum speed, rpm	1300
Battery life	600 hours continuous use
Plate size mm	150 x 150
Overall dimensions (wxdxh), mm	150 x 160 x 70

S512-342



## Stirrers, Magnetic, Big Squid

ika

Compact magnetic stirrer for stirring quantities up to 1.5 litres. Electronically controlled motor allows infinitely adjustable speeds from 0-2500rpm. Excellent chemical resistance due to glass top and synthetic base constructed of TPC-ET.

Catalogue No	<b>S512-342</b>
Model	White
Stirring quantity, litres	1.5
Speed range, rpm	0 - 1500
Plate dimensions, mm	160
Overall dimensions (wxdxh), mm	180 x 203 x 45

S512-710



## Stirrers, Magnetic

stuart

Available in a robust stainless steel or chemically resistant ceramic top. Supplied with fitting for retort stand and 2 x 25mm PTFE coated stirrer bars.

Catalogue No	<b>S512-710</b>	<b>S512-718</b>
Model	CB161	SB161
Plate material	Glass ceramic	Stainless steel
Stirring quantity, litres	15	15
Speed range, rpm	100 - 1500	100 - 1500
Plate dimensions, mm	160 x 160	160 x 160
Overall dimensions (wxdxh), mm	190 x 300 x 110	190 x 300 x 110

S512-734



## Stirrers, Magnetic

Epoxy painted metal case gives excellent resistance to chemicals and corrosion in general. Suitable for the majority of laboratory stirring requirements.

Catalogue No	<b>S512-734</b>	<b>S512-738</b>
Model	AGE	ATE
Speed range	0 - 1200 rpm	0 - 1100 rpm
Max. Stirring Volume	8 litres	25 litres
Power	40W	15W
Voltage	220 - 240v / 50Hz	
Dimensions ( wxdxh), mm	171 x 190 x 75	250 x 285 x 120