

Stoppers

Aa

Bb

Cc

Dd

Ee

Ff

Gg

Hh

Ii

Jj

Kk

Ll

Mm

Nn

Oo

Pp

Qq

Rr

Ss

Tt

Uu

Vv

Ww

Xx

Yy

Zz



Stoppers, Hollow Blown Glass

- Ideal for use with volumetric flasks etc having 'C' length sockets
- Manufactured from Pyrex® borosilicate

Catalogue No	Cone size	Overall height mm
S716-531	10/13	42
S716-532	12/14	47
S716-533	14/15	50
S716-535	19/17	55
S716-536	24/20	61
S716-537	29/22	68



Stoppers, Solid Glass

- Ideally suited for use with Pyrex® bottles and aspirators
- Robust, solid glass construction
- Hexagonal head for easy grip

Catalogue No	Cone size	Head diameter mm	Overall height, mm
S716-572	14/15	32	27
S716-574	19/26	36	38
S716-575	24/29	46	42
S716-576	29/32	56	46
S716-580	45/40	77	55
S716-582	60/46	94	62



Stoppers, Polyethylene

- Ideally suited for use with volumetric flasks and measuring cylinders etc. which have longer 'B' length sockets
- Manufactured from polyethylene
- Large head provides a good grip and gives protection if glassware is accidentally knocked over

Catalogue No	Joint size	Pack qty.
S716-651	10/19	10
S716-652	12/21	10
S716-653	14/23	10
S716-655	19/26	10
S716-656	24/29	10
S716-657	29/32	10
S716-658	34/35	10
S716-659	45/40	10



Stoppers, PTFE

Robust and portable flat head stoppers that fit all standard ground glass jars. The trapezoidal rings provide a high vacuum performance and reduce seizure.

Catalogue No	Joint size	Overall height, mm
S716-682	14/23	34
S716-684	19/26	38
S716-686	24/29	44
S716-688	29/32	47

Stopper, Natural Red Rubber

- Complies with B.S. 2775
- Will withstand most chemical agents in normal laboratory use
- Temperature Range -40 to +70PC

Solid	1-Hole	2-hole	Dimensions (bot x top x l), mm	Pack qty.
S720-003			3 x 4.5 x 12	50
S720-004			4 x 5.8 x 14	50
S720-005	S731-005		5 x 7 x 16	50
S720-006	S731-006		6 x 8 x 16	50
S720-007	S731-007		7 x 9 x 16	20
S720-008	S731-008		8 x 10.5 x 20	20
S720-009	S731-009		9 x 11.5 x 20	20
S720-010	S731-010		10 x 12.5 x 20	20
S720-011	S731-011		11 x 14 x 24	20
S720-013	S731-013		13 x 16 x 24	20
S720-015	S731-015	S742-015	15 x 18 x 24	20
S720-017	S731-017	S742-017	17 x 20.3 x 26	20
S720-019	S731-019	S742-019	19 x 22.5 x 28	20
S720-021	S731-021	S742-021	21 x 24.5 x 28	20
S720-023	S731-023	S742-023	23 x 26.5 x 28	20
S720-025	S731-025	S742-025	25 x 28.5 x 28	20
S720-027	S731-027	S742-027	27 x 31 x 32	10
S720-029	S731-029	S742-029	29 x 33 x 32	10
S720-031	S731-031	S742-031	31 x 36 x 35	10
S720-033	S731-033	S742-033	33 x 38.5 x 38.5	10
S720-035	S731-035	S742-035	35 x 40.5 x 38.5	10
S720-037	S731-037	S742-037	37 x 42.5 x 38.5	10
S720-039	S731-039	S742-039	39 x 44.5 x 38.5	10
S720-041	S731-041	S742-041	41 x 47 x 42	5
S720-043	S731-043	S742-043	43 x 49 x 42	5
S720-045	S731-045	S742-045	45 x 51 x 42	5
S720-049	S731-049	S742-049	49 x 55 x 42	5
S720-053	S731-053	S742-053	53 x 59.5 x 45.5	2
S720-057	S731-057	S742-057	57 x 64 x 49	2
S720-062	S731-062	S742-062	62 x 69 x 49	2
S720-067	S731-067	S742-067	67 x 74 x 49	2
S720-071	S731-071	S742-071	71 x 78 x 49	1
S720-075	S731-075	S742-075	75 x 85 x 54	1
S720-079	S731-079	S742-079	79 x 89 x 54	1
S720-086	S731-086	S742-086	86 x 96 x 54	1
S720-092	S731-092	S742-092	92 x 102 x 54	1
S720-099	S731-099	S742-099	99 x 109 x 54	1
S720-105	S731-105	S742-105	105 x 116.5 x 63.5	1

Stopper, Natural Red Rubber, Solid

Pack of 67 assorted solid red rubber stoppers comprising:

- 5 of each sizes 7, 8, 9, 10, 11, 13, 15, 17, 19, 21, 23
- 2 of each sizes 25, 27, 29, 31, 33, 35

Catalogue No

S743-100

Stopper, Natural Red Rubber, One-Hole

Pack of 67 assorted one-hole red rubber stoppers comprising:

- 5 of each sizes 7, 8, 9, 10, 11, 13, 15, 17, 19, 21, 23
- 2 of each sizes 25, 27, 29, 31, 33, 35

Catalogue No

S743-111

Stopper, Natural Red Rubber, Two-Hole

Pack of 48 assorted two-hole red rubber stoppers comprising:

- 5 of each sizes 15, 17, 19, 21, 23, 25, 27, 29
- 2 of each sizes 31, 33, 35, 37

Catalogue No

S743-222



Aa

Bb

Cc

Dd

Ee

Ff

Gg

Hh

Ii

Jj

Kk

Ll

Mm

Nn

Oo

Pp

Qq

Rr

Ss

Tt

Uu

Vv

Ww

Xx

Yy

Zz

Stoppers

Aa

Bb

Cc

Dd

Ee

Ff

Gg

Hh

Ii

Jj

Kk

Ll

Mm

Nn

Oo

Pp

Qq

Rr

Ss

Tt

Uu

Vv

Ww

Xx

Yy

Zz



Stopper, Silicone to B.S. 2775

- Non toxic, Translucent
- High resistance to many chemicals
- Temperature range -60 to +200PC

Solid	Dimensions (bot x top x l), mm	Pack qty.
S746-003	3 x 4.5 x 12	50
S746-004	4 x 5.8 x 14	50
S746-005	5 x 7 x 16	50
S746-006	6 x 8 x 16	50
S746-007	7 x 9 x 16	50
S746-008	8 x 10.5 x 20	20
S746-009	9 x 11.5 x 20	20
S746-010	10 x 12.5 x 20	20
S746-011	11 x 14 x 24	10
S746-013	13 x 16 x 24	10
S746-015	15 x 18 x 24	10
S746-017	17 x 20.3 x 26	10
S746-019	19 x 22.5 x 28	5
S746-021	21 x 24.5 x 28	5
S746-023	23 x 26.5 x 28	5
S746-025	25 x 28.5 x 28	2
S746-027	27 x 31 x 32	2
S746-029	29 x 33 x 32	2
S746-031	31 x 36 x 35	2
S746-033	33 x 38.5 x 38.5	2
S746-035	35 x 40.5 x 38.5	1
S746-037	37 x 42.5 x 38.5	1
S746-039	39 x 44.5 x 38.5	1
S746-041	41 x 47 x 42	1
S746-043	43 x 49 x 42	1
S746-045	45 x 51 x 42	1
S476-049	49 x 55 x 42	1
S746-053	53 x 59.5 x 45.5	1
S746-057	57 x 64 x 49	1
S746-062	62 x 69 x 49	1
S746-067	67 x 74 x 49	1
S746-071	71 x 78 x 49	1

Stoppers, Rubber, Turnover, Red

With turnover flange and serrations to grip the inside of narrow neck containers. In packs of 10.



Catalogue No	Size, mm	Plug Dia. mm
S746-309	9	6.5
S746-313	13	8
S746-317	17	9.5
S746-321	21	11
S746-325	25	12.5
S746-329	29	14
S746-333	33	16
S746-337	37	17.5
S746-341	41	19
S746-345	45	20.5
S746-349	49	22
S746-353	53	24
S746-357	57	25.5