

Aa

Bb

Cc

Dd

Ee

Ff

Gg

Hh

Ii

Jj

Kk

Ll

Mm

Nn

Oo

Pp

Qq

Rr

Ss

Tt

Uu

Vv

Ww

Xx

Yy

Zz



T352-010

## Clock, Wall Mounting, Quartz Movement

Catalogue No	Case diameter, mm	Colour	Case depth, mm	Dial type
T352-010	255	White	40	White, 12 hr
T352-020	303	White	40	White, 24 hr



T352-060

## Clock, Wall Mounting

A large wall clock featuring temperature and humidity meters. Clear and large display makes it ideal for a variety of settings. Supplied complete with battery.

Catalogue No	Thermometer Range	Humidity Range	Dimensions
T352-060	-35 to 55°C	0 to 100% RH	330 x 280 x 40mm



T352-131

## Timer, Analogue, 1 Hour

Analogue timer counts down from 60 minutes with alarm. Red LED light flashes when counting down. Hanging hook, 60 x 60mm, battery AA.

Catalogue No	T352-131
--------------	----------



T352-165

## Timer, Countdown with Memory

Counts down from 99 minutes, 59 seconds, with alarm. Memory, back clip with magnet, 55 x 55mm, battery AAA.

Catalogue No	T352-165
--------------	----------



T352-172

## Timer, Countup/Countdown with Memory

Counts up and down to/from 24 hours, with alarm. Memory, clock, back clip with magnet, 60 x 50, batteries supplied.

Catalogue No	T352-172
--------------	----------



T352-199

## Timer, Countdown/Countdown with Memory

Counts up and down to/from 10 hours, with 10 key number setting, alarm and clock, magnet on casing. 72 x 65mm, battery AAA.

Catalogue No	Description
T352-199	Timer, countdown/countdown with memory



T352-255

## Timer, Multifunction, Countup

1/100 second stopclock counting up to 10 hours, lap and cumulative split options, three sockets to allow external switching. Black case 100 x 120 x 45mm. Battery R6/AA. Please note that leads are not supplied with this unit.

Catalogue No	T352-255
--------------	----------



T352-265

## Timer, Countdown/Countup, Large Display

Large LCD display, two timing channels counting down 99 hours, 59 minutes, 59 seconds. Count-up facility in normal time. Start, stop and re-set functions. White case 100 x 120 x 46mm, R6/AA battery.

Catalogue No	T352-265
--------------	----------



T352-270

## Timer, Countdown/Countup, Large Display, Decimal Seconds

Large LCD display counting up with 4 timing channels; seconds, minutes, hours and decimal seconds. Count down timer with memory and alarm at the end of count down period. Normal time with calendar and alarm. Case 90 x 160 x 70mm. R14/C size battery.

Catalogue No	T352-270
--------------	----------

## Timer, Countdown/Countdown with Memory

Single channel bench timer, 11mm LCD display. Countdown mode with alarm on expiry. Count-up mode until stopped. Memory for countdown times. Case 85 x 85 x 45mm. Supplied with battery.

Catalogue No **T352-282**

T352-282



## Timer, Countdown/Countup, Large Display

16 jewel lever movement, 1/5th second, 60 minute, 60 second, chrome finish metal case. 52 mm diameter, stop/start top crown button, side hand re-set.

Catalogue No **T352-412**

T352-412



## Timer, Countdown/Countup, Large Display,

1/100 second stopwatch with lap facility, time and calendar, alarm and hourly chime, 9mm LCD display, water resistant cases 66 x 52mm, with neck cord. Batteries included.

Catalogue No **T352-430**

T352-430



## Decimal Seconds

1/100 second stopwatch with lap facility, time and calendar with 12/24 hour display and alarm. 6mm digit height, water resistant case 70 x 62mm with neck cord. Battery CR2032.

Catalogue No **T352-444**

T352-444



## Timer, Bench Process

1/100 second stopwatch with lap and cumulative lap and split facility, water resistant case 71 x 63mm, with neck cord. Batteries included.

Catalogue No **T352-456**

T352-456



## Note:

For Talley Counters please refer to **page 79**



Aa

Bb

Cc

Dd

Ee

Ff

Gg

Hh

Ii

Jj

Kk

Ll

Mm

Nn

Oo

Pp

Qq

Rr

Ss

**Tt**

Uu

Vv

Ww

Xx

Yy

Zz