

Tongs, Trays

Aa
Bb
Cc
Dd
Ee
Ff
Gg
Hh
Ii
Jj
Kk
Ll
Mm
Nn
Oo
Pp
Qq
Rr
Ss
Tt
Uu
Vv
Ww
Xx
Yy
Zz



T553-014



T553-018

Tongs, Crucible, Length 200mm

Catalogue No	Material	Pattern
T553-012	Mild steel	Straight
T553-014	Stainless steel	Straight
T553-016	Mild steel	Bow
T553-018	Stainless steel	Bow



T609-264

Trays, Instrument, Polypropylene

- Rigid and robust
- Autoclavable
- Vertical sides

Catalogue No	Dimensions (lxwxh), mm
T609-260	350 x 200 x 22
T609-264	350 x 250 x 22
T609-268	610 x 350 x 22



T609-285

Trays, Laboratory, White, PVC

- Lightweight and unbreakable
- End carrying handles

Catalogue No	Dimensions (lxwxh), mm
T609-281	422 x 321 x 24
T609-283	410 x 321 x 25
T609-285	524 x 236 x 27
T609-287	555 x 406 x 24



T609-320

Trays, Tidy, Polystyrene

- Lightweight, compartmentalised trays
- All trays are 340 x 440mm
- Designed primarily to fit into and remain inside a drawer

Catalogue No	No of compartments	Compartment size (lxwxh), mm
T609-318	3	2 x (100 x 400 x 55 deep) 1 x (140 x 400 x 55 deep)
T609-320	12	12 x (95 x 95 x 55 deep)
T609-324	5	2 x (90 x 190 x 65 deep) 1 x (95 x 400 x 65 deep) 2 x (120 x 190 x 65 deep)



T609-350

Trays, White, PVC

Sloping sides enable the dishes to be stacked. PVC can be used at 75°C for short periods, it is however recommended that 65°C is not exceeded for longer periods.

Catalogue No	Dimensions
T609-350	150 x 200 x 45
T609-352	200 x 250 x 60
T609-354	260 x 320 x 70
T609-356	300 x 350 x 80
T609-358	330 x 430 x 90
T609-360	420 x 520 x 90
T609-362	530 x 675 x 100
T609-364	420 x 540 x 180