



Wagtech Court
Station Road
Thatcham
Berkshire RG19 4HZ
United Kingdom

Tel +44 (0) 1635 872929
Fax +44 (0) 1635 862898
info@wagtech.co.uk
www.wagtech.co.uk

Landrover 110

Mobile Laboratories Specifications



Wagtech Mobile Laboratory

Land Rover Defender Left Hand Drive 110 Station Wagon, 2402cc Diesel engine

VEHICLE SPECS:

Export tropical specification

- 2.4 litre diesel engine (111 bhp)
- EU2 Emissions - EU4 also available depending on market requirement
- Hot climate cooling - oil cooler
- Raised Air Intake
- Power assisted steering
- 6 speed manual gearbox with locking differential
- Permanent 4-wheel drive
- 5 door with 5 seats in Vinyl trim, cab and second row seats retained
- Front and rear 3 point seat belts
- Colour - White - several other colours available if required. Metallic paint carries an additional cost.
- Payload 1183 kg, gww 3500 kg
- Aluminium corrosion resistant body
- Laminated windscreen
- Coil spring suspension
- Heavy duty double acting shock absorbers
- Disc brakes all-round
- Michelin 4 x 4 - 7.50x16 tubed radial on/off road tyres
- Spare wheel assembly on rear door
- Locking filler cap
- Rear mud flaps
- Exterior rear view mirrors
- 16.5 gall (75 litre) fuel tank
- 2 speed wipers with intermittent wipe
- Heater/demister & rear wash wipe
- Underside protection bar
- Rubber mats - front foot well and second row
- Tool kit
- FM/MW radio/single play CD
- Roof rack with access ladder
- Air-conditioning
- Front nudge bar (Grill guard) - CODAN type with aerial mounting brackets
- Pre-delivery inspection - carried out on all vehicles
- Supplied with a 36 month or 100,000km International Warranty - whichever comes first
- Sign written to requirements

Rear Fit-Out Specification:

Specially modified to suit the purpose of a general laboratory. Consists of the following factory fitted components:

- Exterior Grade ply floor overlaid with washable anti-slip flooring
- Sink top and waste with clean and dirty water containers to left hand side of rear saloon
- 12V Fridge/Freezer: Capacity 43litres including 4 litre freezer compartment under laboratory laminated work surface at rear of second row seats. Powered by 12V vehicle electronics.
- Lockable cupboard with shelf left of fridge.
- Front opening cupboard to right of refrigerator over wheelbox.
- Separate intermediate battery with split charge relay.
- Hook up for generator (generator not provided)
- Rear overhead lighting.



Land Rover 110

Illustrations showing alternative configuration on standard 110 model

